

Minor Surgery Tray

Catalog # 3555/ *

Skin Biopsy Kit

Disposable tray and instruments eliminate risk of cross-infection.

Specially designed for emergency and physician office use.

Contains extensive selection of items needed for surgical procedures.

Can be used for most types of emergency and small surgical procedures.

Designed for removal of foreign bodies, or occlusion of bleeding vessels.

Contains extensive selection of items required to perform skin biopsy procedure.

Sutures can be incorporated.



Standard Tray (# 3555)

Europlas® Fine point scissors 4¾"

Europlas® Hemostat curved 5½"

Europlas® Needle holder 5½"

Europlas® Sharp/blunt scissors 4¾"

Europlas® Lift Forceps 5"

Europlas® Forceps deluxe forceps fine point 4"

1 pair Latex gloves

Fenestrated drape

Polybacked towel

4 gauze sponges 10X10 cm (4"x4")

5 round gauze sponges

2 Sutures 3:0

2 scalpels #. 11 & 15

Stainless steel forceps 4"

Alcohol prep pad

Povidone iodine prep pad

Waste bag with adhesive tape



* Other Compositions according to customer request

Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

