

Circumcision Tray

Catalog # 931 & 932

Convenient, Disposable, Sterile

Contains all the necessary components for complete procedure.

Clinically proven **Eurofor**[®] Mogen Circumcision Clamp.

Mogen Circumcision Clamp prevents side slipping protecting Glans Penis.

Mogen Clamp ensures accurate visibility and gauging of Prepuce and Mucosa to be excised.

Europlas[®] Probe for easy separation of Prepuce from Glans Penis.

Disposable tray and instruments eliminate risk of cross infection.



Complete Circumcision Tray (#931)

Eurofor[®] Mogen circumcision clamp **Europlas**[®] Probe

Stainless steel hemostat

Europlas[®] Forceps

Europlas[®] Sharp/blunt scissors

Disposable scalpel

3 gauze sponges

Paraffin tape

1 pair latex gloves

Povidone iodine swab

Cotton tipped applicator

Waste bag with adhesive tape

Poly backed towel

Fenestrated drape



Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.

Tel: 972 3 9022885

Fax: 972 3 9024531

euroband@hotmail.com

www.euroband.com

EUROBAND

