

CYSTOSCOPY PREP TRAY

A Sterile Tray Ideally Designed to Prepare the Patient for Cystoscopy Procedure.

Rigid PVC transparent tray; Built-in separate compartment for liquids; Peel-away Tyvek® cover.

Components:

- Absorbent poly-backed drape for use as sterile field
- Fenestrated drape for use around body member
- 1 pair latex gloves
- 3 gauze balls for use with disinfectant
- 1 lift forceps
- 3 (10x10cm) gauze sponges
- 1 sachet antiseptic solution
- PENICLAMP™ - an effective device - holds the male meatus in "closed" position, preventing leakage after lubricant has been introduced in preparation for insertion of cystoscope.



- Syringe, ** pre-filled with Pedicat™ lubricant.
**The 6ml syringe features:
 - a) Tinted airtight transparent syringe protects contents from light.
 - b) Specially designed protective cap preserves sterility up until the time of use;
 - c) Innovative rotating piston handle. Prior to use, rotate the piston handle as indicated, to ensure smooth, soft, and easy release of the piston.



CYSTOSCOPY
PREP TRAY

Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

