

Diapouch™


WATERPROOF POUCH THAT HOLDS AND PRESERVES CATHETER TUBES

Sterile, Elastic & Soft Waterproof Pouch
Outer size: 17X8cm, inner size: 13X7cm.

Holds and preserves catheter tubes and connectors of
Hemodialysis/Peritoneal/Oncology catheter

During showering, provides fair protection against water penetration into
the pouch and around the exit site*

Protects and manages entry/exit sites of the catheter and its tubing and
connectors.



CATHETER IN DIAPOUCH™



CATHETER PRESERVATION WATERPROOF PACK*

Standard Products:

DIAPOUCH™ - Catalog No. 01140:

Sterile, packed in flexible transparent peel-open pouch

CATHETER PRESERVATION PACK - Catalog No. 01143:

Sterile, packed in flexible transparent peel-open pouch

Contains:

- **DIAPOUCH™** - Waterproof Pouch outer size: 17 X 8 cm
- **ADHESIVE PLASTER** air permeable, size:10X12 cm

***Note:** *The pouch is completely water resistant. However, pouch flap should completely cover the absorbent dressing around the exit site, with no extruding edges. The adhesive tape should be applied flatly with no kinks and cover the entire flap.*

Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.

Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com

www.euroband.com

EUROBAND

