

MICRO100™ PROTECTIVE HOOD MASK

- ◆ Electrostatic filter blocks 99.9999 % of microorganisms and viruses, and 99.7 % of particles
- ◆ Covers all breathing orifices, and protects ears and eyes
- ◆ Seal around neck prevents penetration of contaminated air, allowing respiration via filter only
- ◆ Additional layer contains antiseptic properties - microorganisms in the filter are destroyed
- ◆ The Micro100™ filter is designed to provide excellent breathing performance, until clogged by heavy dust or particles
- ◆ Tests reveal that when Micro100™ is properly worn, with its filter opposite the mouth, there is limited Carbon Dioxide (CO2) accumulation within the hood
- ◆ One size fits all
- ◆ Hood is easily removed. The “bag shape” is removed by pulling from the bottom upwards - to be turned inside-out to accommodate also discarded disposable garments and gloves
- ◆ After removal, when contaminated surfaces are unexposed, the hood mask can be discarded safely

To safeguard workers in contaminated areas or vulnerable environments such as chicken huts, maximum protection of respiratory orifices and eyes is advised.

Micro100™ is a disposable device providing the optimal solution to ensure ultimate security at minimum expense.

Disposable - Maximum Protection - Minimum Expense.

POLLAK (INTERNATIONAL) LTD.- 20 Hamlacha St. Afek, Rosh Ha'ayin 48091 Israel. www.euroband.com.

Tel:+972 3 902 2885 Fax:+972 3 902 4531. info@muski.com

Complies to ISO 9001:2000; ISO 13485:2003; FDA registered facility & CE marked products.

Pat. Pending