

Peritoneal/Oncology Dressing Set


Catalog # 1144

Sterile waterproof dressing* & pouch, protects exit site of peritoneal/oncology catheter; holds and preserves catheter tubes and connectors.

Includes:

- **Silvalpad™** – Absorbent antibacterial dressing.
- **Diapouch™** - Elastic waterproof Pouch for Catheter Connectors. Outer-size; 17 cm. X 8 cm.
- **Adhesive tape** 10 cm. X 12 cm.

*** Note:** The pouch is completely water resistant. However, pouch flap should completely cover the absorbent dressing around the exit site, with no extruding edges. The adhesive tape should be applied flatly with no kinks and cover the entire flap.



Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

