

# Tracheostomy Care Tray

Catalog # 921/★

Tray Designed  
for Convenience  
in Cleaning and  
Care of  
Tracheostomized  
Patients.



## Standard Kit Components:

- ❖ Large tray, with 3 division extra deep compartments for holding liquids and working components.
  - ❖ Drape - serves as a sterile field.
  - ❖ Comfortable, form-fitting latex gloves.
  - ❖ Stainless steel & nylon<sup>®</sup> brush prevents corrosion.
  - ❖ Smooth-tip design guards against potential damage to tracheostomy tube.
  - ❖ 75 cm Medical-grade twill tape contains no irritants.
  - ❖ Tracheostomy 10x10 cm Non-Adherent dressing.
  - ❖ 2 Plastic, cotton-tipped applicators.
  - ❖ 2 Large pipe cleaners.
  - ❖ 4 Gauze sponges 10x10cm.
  - ❖ Optional: Lift forceps for sponges to assist in cleaning procedure.
  - ❖ Optional: Pair of latex gloves.
- ★ Available in different versions.



**Pollak Ltd.** Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.  
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com  
[www.euroband.com](http://www.euroband.com)

**EUROBAND**

