

Pedicat™

Catalog # 881

Catheter Lubricating Gel

**WATER SOLUBLE
ANTISEPTIC
LUBRICATING GEL.**


**Aluminum Laminated
Sachet -**

**The sachet can be opened
in two ways –**

- Peel apart at top edge - catheter tip can be inserted into the sachet for lubrication
- Tear open at bottom corner, to enable the gel to be squeezed out on/into the organ for lubrication.



When tearing open at bottom corner, the gel can be squeezed out on/into the organ for lubrication.



Peel apart at top edge – catheter tip is inserted into the sachet for lubrication.

Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

